

Planarian Scrunching as a Quantitative Behavioral Readout for Noxious Stimuli Sensing

Ziad Sabry¹, Christina Rabeler¹, Danielle Ireland¹, Kevin Bayingana¹, Eva-Maria S. Collins^{1,2}

¹ Biology Department, Swarthmore College ² Physics Department, UC San Diego

Corresponding Author

Eva-Maria S. Collins
ecollin3@swarthmore.edu

Citation

Sabry, Z., Rabeler, C., Ireland, D., Bayingana, K., Collins, E.M.S. Planarian Scrunching as a Quantitative Behavioral Readout for Noxious Stimuli Sensing. *J. Vis. Exp.* (161), e61549, doi:10.3791/61549 (2020).

Date Published

July 30, 2020

DOI

10.3791/61549

URL

jove.com/video/61549

Abstract

Freshwater planarians normally glide smoothly through ciliary propulsion on their ventral side. Certain environmental conditions, however, can induce musculature-driven forms of locomotion: peristalsis or scrunching. While peristalsis results from a ciliary defect, scrunching is independent of cilia function and is a specific response to certain stimuli, including amputation, noxious temperature, extreme pH, and ethanol. Thus, these two musculature-driven gaits are mechanistically distinct. However, they can be difficult to distinguish qualitatively. Here, we provide a protocol for inducing scrunching using various physical and chemical stimuli. We detail the quantitative characterization of scrunching, which can be used to distinguish it from peristalsis and gliding, using freely available software. Since scrunching is a universal planarian gait, albeit with characteristic species-specific differences, this protocol can be broadly applied to all species of planarians, when using appropriate considerations. To demonstrate this, we compare the response of the two most popular planarian species used in behavioral research, *Dugesia japonica* and *Schmidtea mediterranea*, to the same set of physical and chemical stimuli. Furthermore, the specificity of scrunching allows this protocol to be used in conjunction with RNA interference and/or pharmacological exposure to dissect the molecular targets and neuronal circuits involved, potentially providing mechanistic insight into important aspects of nociception and neuromuscular communication.

Introduction

In addition to their popularity for stem cell and regeneration research^{1,2,3}, freshwater planarians have long been used in behavioral studies^{4,5}, taking advantage of their comparatively large size (a few millimeters in length), ease and low cost of laboratory maintenance, and broad

spectrum of observable behaviors. The introduction of computer vision and automated tracking to planarian behavior studies^{6,7,8,9,10,11} have enabled quantitative differentiation of behavioral phenotypes. Animal behavior is a direct readout of neuronal function. Because the planarian

nervous system is of medium size and complexity, but shares conserved key elements with the vertebrate brain^{12, 13, 14}, studying planarian behavior can provide insight into conserved mechanisms of neuronal action which may be hard to directly probe in more complex organisms. Thus, planarians are a valuable model for comparative neurobiology studies^{8, 12, 15, 16, 17, 18, 19, 20, 21}. In addition, the aquatic environment allows for rapid and facile exposure to chemicals to study their effect on brain function in regenerating and adult planarians, making them a popular system for neurotoxicology^{22, 23, 24, 25, 26}.

Planarians possess three distinct gaits, referred to as gliding, peristalsis, and scrunching. Each gait is exhibited under specific circumstances: gliding is the default gait, peristalsis occurs when ciliary function is compromised^{7, 27}, and scrunching is an escape gait – independent of cilia function – in response to certain noxious stimuli⁷. We have shown that scrunching is a specific response, elicited by the sensation of certain chemical or physical cues, including extreme temperatures or pH, mechanical injury, or specific chemical inducers, and thus is not a general stress response^{7, 28, 29}.

Because of its specificity and stereotypical parameters, which can easily be quantified using this protocol, scrunching is a powerful behavioral phenotype that enables researchers to perform mechanistic studies dissecting sensory pathways and neuronal control of behavior^{25, 28}. Additionally, scrunching has been shown to be a sensitive endpoint to assay adverse chemical effects on nervous system development and function in neurotoxicology studies^{22, 24, 25, 30}. As several different sensory pathways seem to converge to induce scrunching through various mechanisms²⁸, scrunching differs from other planarian behaviors because various, but specific, stimuli can be used

to dissect distinct neuronal circuits and study how different signals are integrated to produce the scrunching phenotype.

Importantly, species differences exist, wherein one chemical may elicit scrunching in one planarian species, but a different behavioral response in another. For example, we have found that anandamide induces scrunching in the planarian species *Dugesia japonica* but induces peristalsis in *Schmidtea mediterranea*²⁸. This example highlights the importance of being able to reliably distinguish between the different gaits, because they are the phenotypic manifestations of distinct molecular mechanisms. However, distinction of scrunching from peristalsis is difficult using qualitative observational data, because both gaits are musculature-driven and share qualitative similarities^{7, 28}. Thus, to distinguish the gaits it is necessary to perform cilia imaging or a quantitative behavioral study, which allows distinction based on characteristic parameters^{7, 28}. Because cilia imaging is experimentally challenging and requires specialized equipment such as a high-magnification compound microscope and a high-speed camera^{7, 28}, it is not as broadly accessible to researchers as quantitative behavioral analysis.

Here, we present a protocol for (1) the induction of scrunching using various physical (noxious temperature, amputation, near-UV light) and chemical (allyl isothiocyanate (AITC), cinnamaldehyde) stimuli and (2) the quantitative analysis of planarian behavior using freely available software. By quantifying four parameters (frequency of body length oscillations, relative speed, maximum amplitude, and asymmetry of body elongation and contraction)⁷, scrunching can be differentiated from gliding, peristalsis, and other behavioral states reported in the literature, such as snake-like locomotion¹⁵ or epilepsies¹⁵. Furthermore, while scrunching is conserved among different planarian species⁷, each

species has its own characteristic frequency and speed; therefore, once the gliding and scrunching speeds of a species have been determined, speed alone can be used as a means to distinguish scrunching from gliding and peristalsis²⁹. The protocol assumes no prior training in computational image analysis or behavioral studies and thus can also be applied for planarian behavioral experiments in a teaching laboratory context at the undergraduate level. Example data to facilitate protocol adaptation is provided in the Supplemental Material.

Protocol

1. Quantitative planarian behavior assays

1. Experimental setup

1. Place a dimmable LED panel upon a flat surface.

The LED panel serves two purposes: (1) to provide a uniform white background and (2) to be used as an adjustable light source to obtain appropriate contrast. Place a 100 mm Petri dish arena upon the LED panel.

NOTE: To increase throughput, a multi-well plate may be used as an arena^{23, 24}, but larger arenas facilitate automated image analysis.

2. Mount a camera on a ring stand above the arena (**Figure 1A**). Adjust the camera position, height, and focus as necessary so that the entire arena is centered within the field of view and is in focus (**Figure 1B**).

NOTE: The camera resolution needs to be high enough to clearly distinguish a planarian from the homogenous background provided by the LED panel.

3. Fill the arena with the appropriate exposure media (planarian water or chemical solution) to half-

maximum volume (this will be referred to as a bath). This corresponds to approximately 25 mL for a 100 mm Petri dish. Turn on the LED panel and turn off any other light sources that may negatively affect recording quality (i.e., nearby light sources that produce a glare onto the arena).

CAUTION: Manage hazardous chemical solutions appropriately by wearing full personal protective equipment (PPE) and moving the experimental setup to a fume hood if necessary. Follow federal and state regulations on waste disposal.

4. Drop a planarian toward the center of the arena using a transfer pipette. Begin recording. Record data as image sequences in a native Fiji³¹ format (TIFF, GIF, JPEG, PNG, DICOM, BMP, PGM, or FITS; see image analysis section 1.2).

NOTE: Because behaviors and sensitivity to external stimuli vary among individual planarians, it is important to collect data on a sufficiently large number of biological replicates, in addition to performing technical replicates. We have worked with up to 10 medium-sized (4-7 mm) planarians in a 100 mm Petri dish at once. While time efficient, multiple planarians in the Petri dish at once make data analysis more difficult since planarians may cross paths.

1. For gliding experiments, record using at least 1 frame per second (FPS). For scrunching/peristalsis experiments, record using an FPS that is at least twice the scrunching/peristalsis frequency of the planarian species. If the planarian species has an unknown scrunching/peristalsis frequency, use 10 FPS as a starting point and increase/decrease as appropriate.

2. When using a chemical solution, transfer the planarian using as few drops of planarian water as possible so that the concentration of the chemical solution is not significantly changed.

5. For gliding experiments, record 1-2 minutes of gliding behavior. For scrunching/peristalsis experiments, record long enough to capture at least 3 consecutive oscillations occurring in a straight line. Once the experiment is completed, terminate the recording.

NOTE: For scrunching/peristalsis experiments, if a planarian does not satisfy the termination criterion within a fixed time period that needs to be consistent across replicates and is empirically determined based on the stimulus, terminate the recording and test another planarian.

1. If the planarian reaches the boundary of the arena without satisfying the termination criterion, pipette the planarian back to the center of the arena.

NOTE: Avoid repeated pipetting of an individual for recording, as this may change its behavior.

6. Remove the planarian(s) from the arena and dispose of the planarian water or chemical solution in appropriate waste containers. Planarians that were in planarian water can be returned to their home container.

NOTE: Avoid cross contamination by using different arenas for different media (i.e., gliding in planarian water experiments should not be run in an arena previously used for scrunching/peristalsis experiments with chemical exposure).

1. Serially rinse planarians exposed to a chemical solution in 3 clean 100 mm Petri dishes filled with 25 mL of planarian water to thoroughly dilute out

any chemicals. If scrunching or peristalsis was induced, place these planarians in a separate container. Planarians can be returned to their home container after one month since most cells would have turned over by that time¹.

NOTE: If multiple different experiments are needed for the same population of planarians, e.g., for an RNAi population, allow planarians to recover for 24 hours before running the next experiment. Order the experiments such that the least invasive experiment is first and the most invasive experiment (e.g., amputation) is run last.

2. If running multiple experiments in the same arena, properly dispose of the bath solution and remove any mucus trails by wiping down the arena with a paper towel between runs.

NOTE: The protocol can be paused here.

2. Quantitative analysis of planarian behavior

1. Perform planarian behavior assays as described in Section 1.1.

2. Open the raw image sequence for an experiment in Fiji (**File > Import > Image Sequence**) and select the first image in the image sequence. In the Sequence Options window, check the box for "Sort names numerically" and click "OK". Once the image sequence is loaded in, convert the image sequence to 8-bit (**Image > Type > 8-bit**) and use the arrow tool or slider at the bottom of the image stack to watch or pan through the image sequence.

NOTE: For gliding experiments, all data can be used as long as the planarian can be clearly seen throughout the recording. However, it is usually sufficient to analyze the free motion in the center

of the arena by extracting the relevant part(s) as described below.

3. To extract a time period and region of interest, draw a region of interest encompassing the full path of a planarian using the rectangle tool (**Figure 2A, 2B**). Right click on the image stack and select **Duplicate...**, check the box for **Duplicate stack**, enter the first and last frames of the sequence of interest, and click **OK**. If multiple planarians were imaged simultaneously, repeat this region selection and duplication step for each planarian in the arena so that there are as many open image stacks as there are planarians in the arena. The following steps (Steps 1.2.4-1.2.10) should be performed on each image stack, one at a time.

1. For gliding experiments, extract a period of gliding where the planarian moves at least twice its body length.

NOTE: The more gliding data extracted per planarian, the more reliable the data will be. The planarian does not need to be moving in a straight line for the gliding analysis.

2. For scrunching/peristalsis experiments, extract an instance when the planarian undergoes a minimum of three consecutive (ideally more) body oscillations in a straight line, making sure each oscillation is a complete elongation-contraction cycle, as full oscillations are necessary to accurately determine the frequency.

NOTE: The more oscillations that can be extracted, the more reliable the data will be. Do not use sequences where the planarian is turning as these will result in inaccurate length measurements.

4. Apply a threshold to the duplicated image stack (**Image > Adjust > Threshold**) to binarize the image and extract the planarian from the background. Adjust the sliding bars as necessary such that the entire planarian is highlighted in red. The exact values are dependent on imaging quality. Leave the boxes for **Dark background**, **Stack histogram**, and **Don't reset range** unchecked. Scroll through the image stack to ensure a good threshold range (i.e., the planarian is well separated from the background throughout the stack), and then click **Apply**.

5. In the **Convert Stack to Binary** window, set the Method to **Default** and the Background to **Light**. Uncheck all boxes in this window and then click **OK**. A binarized image showing a black planarian on a white background will appear (**Figure 2C**). Make sure that the entire planarian is visible in all frames of the image sequence.

NOTE: Unwanted objects in the binarized image sequence that are smaller or larger than the planarian can be filtered out in the subsequent analysis using a size filter (**Figure 2Ciii**).

6. Set measurements by clicking **Analyze > Set Measurements**. Check the boxes for **Area**, **Center of mass**, **Stack position**, and **Fit ellipse** and click **OK**.

NOTE: These parameters only need to be set once per Fiji session.

7. Select the open image stack and select **Analyze > Analyze Particles**.

8. In the **Analyze Particles** window, select **Show > Masks** to open a new stack showing all the objects that were detected with the chosen parameters. This can be used to visually check that only measurements

of the planarian are being taken. A size filter may be set at this step to remove unwanted noise by entering the approximate area of the planarian (in pixel² units) in the space provided. Check the boxes for **Display results** and **Clear results** and click **OK**.

NOTE: In the **Results** window, if the index (first column) does not equal the slice number for all rows, this means that either too many or too few objects were tracked. One possibility for this discrepancy is the presence of other objects besides the planarian or that the planarian was not tracked in specific frames.

9. Pan through the mask image stack using the slider at the bottom of the panel. If there is any noise or there are frames that lack a planarian, close the **Results** window and the mask image stack. Repeat steps 1.2.7-1.2.8 by adjusting the area filter to only remove objects other than the planarian.

NOTE: If the planarian is missing from the frame in the mask, this suggests that the lower bound of the area filter was set too high.

10. On the **Results** window, save the data using **File>Save As**. Add the .csv extension to the filename to save data as comma-separated values. Once data for the image stack is saved, close the respective image stack, and **Results** and **Mask** windows.

11. Import data and further analyze using any spreadsheet software or freeware. To calculate gliding speed, refer to section 1.3. To calculate the scrunching/peristalsis full parameter set, refer to section 1.4.

NOTE: The protocol can be paused here.

12. To determine the pixel to actual length conversion, open an image in Fiji with a reference length (e.g., the

diameter of the arena). Select the line tool and draw a line over the known length.

13. Convert pixel units to a standard unit of length by clicking **Analyze > Set Scale**. Enter the length corresponding to the line drawn on the image in the **Known distance** box and change **Unit of length** from **pixel** to the chosen standard unit of length. The conversion factor is written next to **Scale**.

NOTE: A pixel conversion value is not required for gliding or scrunching/peristalsis analyses in sections 1.3 and 1.4.

3. Calculation of gliding speed

1. Using the data file saved in Section 1.2, load the center of mass (COM) x and y coordinates and the major axis data. If the data is saved as a comma-separated values file, these lists correspond to the "XM", "YM", and "Major" columns, respectively.

2. Calculate the displacement (*d*) of the planarian center of mass in pixels for each frame with respect to the next frame using the "XM" and "YM" data columns. Displacement (*d*) is given by:

$$d = \sqrt{(x_2 - x_1)^2 + (y_2 - y_1)^2}$$

where *x*₁ and *y*₁ refer to the COM coordinates (XM, YM) of one frame and *x*₂ and *y*₂ refer to the COM coordinates (XM, YM) of the subsequent frame.

3. Set the planarian body length as the 95th percentile of the "Major" column. Since planarians exhibit a wall preference behavior³², this ensures that the calculated planarian body length is representative of when the planarian is elongated²⁴.
4. Normalize displacement by planarian body length by dividing the pixel displacements per frame by the

planarian body length (l). Normalized displacement (d_n) is given by:

$$d_n = \frac{d}{l}$$

5. Generate a list of normalized speeds by dividing the normalized displacements by the time elapsed per frame (inverse of the recorded FPS). Normalized gliding speed (s_n) is given by:

$$s_n = \frac{d_n}{(FPS)^{-1}}$$

6. Calculate the normalized gliding speed of the planarian by taking the average of the normalized speeds list (s_n). The standard deviation may be used as an uncertainty measurement for the planarian.
7. Repeat steps 1.3.1-1.3.6 for each planarian to be analyzed. Average and take the standard deviation of the gliding speeds for all planarians to get the gliding speed and associated uncertainty, respectively, for a planarian population.

4. Distinction of scrunching and peristalsis gaits using the full parameter set

1. Load the major axis data list from the data file saved from Section 1.2. If the data is saved as a comma-separated values file, this corresponds to the **Major** column.
2. Create a list that numbers each data point in the **Major** column, starting with 0. Convert this list to time elapsed per frame by dividing by the recorded FPS.
3. Plot the **Major** column data with respect to time elapsed to generate a scrunching/peristalsis oscillation plot (**Figure 3A**). Using the oscillation plot, trim the data to at least three consecutive, straight-line oscillations (**Figure 3Bi**). Trim the data

to start and end at local peaks (maximum elongation of oscillation) or troughs (minimum elongation of oscillation).

NOTE: If local extrema are not approximately equal (peaks/troughs differ dramatically in heights), this suggests that the oscillations are not straight-line (**Figure 3Bii**). Extract another sequence of at least three consecutive, straight-line oscillations. Refer to Section 1.2.

4. Confirm that the oscillation sequence of interest has been extracted and trimmed properly by replotting the trimmed **Major** data with respect to time. Use this trimmed data list for all subsequent calculations.
5. To calculate oscillation frequency (v_m), divide the number of oscillations (O_n) by the total number of data points in the trimmed major axis data list (N). Multiply FPS by this value to get frequency in oscillations per second.

$$v_m = FPS * \frac{O_n}{N}$$

6. To calculate maximum elongation ($|\Delta\epsilon|_{max}$), subtract the absolute minimum body length (l_{min}) from the absolute maximum body length (l_{max}). Normalize to elongated body length by dividing by the absolute maximum body length.

$$|\Delta\epsilon|_{max} = \frac{l_{max} - l_{min}}{l_{max}}$$

7. To calculate speed per body length (v_m^*), multiply the calculated maximum elongation by the oscillation frequency.

$$v_m^* = v_m * |\Delta\epsilon|_{max}$$

NOTE: Speed alone can be used to distinguish between scrunching and peristalsis gaits⁷.

8. To calculate the fraction of time spent elongating (f_{elong}), take the derivative of the trimmed major axis data list with respect to time. Divide the number of positive data points (i.e., when the derivative is >0 (n_p), by the total number of data points in the major axis data list (n_t)).

$$f_{elong} = \frac{n_p}{n_t}$$

NOTE: Scrunching planarians exhibit an asymmetric fraction of time spent elongating whereas planarians performing peristalsis spend equal amounts of time elongating and contracting⁷.

9. Repeat steps 1.4.1-1.4.8 for each planarian to be analyzed. Calculate a planarian population parameter set by taking the average and standard deviation of each parameter.

NOTE: The parameter set can be used to determine if the oscillation behavior is scrunching, peristalsis or some other form of locomotion with periodic body shape changes. Both scrunching and peristalsis have fixed parameters for a given species⁷, with scrunching parameters generally being greater than peristalsis parameters⁷. While it is possible that one of the parameters may fall outside of the species-specific range, as we have previously observed with chemical induction²⁸, the observed behavior must agree with at least 3 of 4 published parameters to be categorized as either peristalsis or scrunching.

2. Scrunching induction

1. Physical stimuli (noxious temperature, UV light, amputation)

1. For all physical stimuli experiments, refer to Section 1.1 for the experimental setup.

NOTE: It is best to use a large arena, such as a 100 mm Petri dish, for physical stimuli experiments to allow for more open space for maneuvering a pipette and/or razor blade.

2. To induce scrunching via noxious temperature, heat planarian water in a glass beaker (at least 100 μ L per planarian to be tested) to 65 °C on a hot plate.

1. Place a planarian in the center of the arena. Wait until the planarian orients itself upright and begins gliding. Begin recording.

2. Using a P-200 pipette, slowly pipette 100 μ L of the 65 °C planarian water post-pharyngeally onto the tail end of the planarian to induce scrunching.

NOTE: Make sure the heated planarian water stays at 65°C. If necessary, reheat the water to 65°C prior to starting another experiment. Since pressure can also induce scrunching, slow pipetting is necessary. Pipetting room temperature water in the same way as in the experiment can serve as a control and practice option.

3. Stop the recording once scrunching has ceased. Place the planarian in a recovery container and exchange the media in the petri dish with fresh, room temperature planarian water if running more experiments.

3. To induce scrunching via amputation, transfer a planarian to the center of the arena and wait until the planarian orients itself upright and begins gliding. Begin recording.

1. Amputate the planarian using a clean razor blade. Amputations may be done anywhere along the planarian as long as the cut location is consistent across experiments.

NOTE: Scrunching parameters are extracted from the anterior piece. Thus, avoid obstructing the camera's view of this part of the planarian when applying the cut by approaching from the posterior end. Plastic cover slips also work well for cutting and are a safer option, especially in a teaching setting.

2. Stop the recording once the anterior piece has ceased scrunching. Remove both pieces, place them in a separate container and allow them to regenerate for 7 days. Amputated planarians can be reincorporated into the home container once regenerated.

4. To induce scrunching using near-UV light, attach appropriate filters (e.g., Roscolux filters) to the camera lens to reduce the amount of reflected near-UV light that is collected by the camera and may interfere with imaging the planarian's response. Instead of using the LED panel to illuminate the arena from below, use ambient red lighting to which planarians are insensitive³³.

1. Fill a 100 mm Petri dish arena with planarian water and place a single planarian (5-9 mm) in the center of the arena. Begin recording at 10 FPS.
2. Hold a Class II UV laser pointer (405 ± 10 nm, output power <5 mW) approximately 30 cm from the arena. Position the laser pointer at a 45° angle from the gliding planarian and then shine the laser pointer for 5-10 seconds halfway between the

posterior end of the pharynx and the tail tip to induce scrunching.

NOTE: The power of the laser pointer can be measured using a near-UV-sensitive power meter.

3. Wait for the planarian to start gliding again before attempting two more stimulations on the same individual to test for reproducibility of the reaction. If the planarian keeps showing the same behavior, stop recording and put the planarian back in its container. If the behavior changes between stimulations, additional tests will show which response is the most prominent.

NOTE: Planarians can become desensitized to near-UV light and will stop reacting. Consecutive stimulations require a rest period of 8-10 seconds.

2. Chemical stimulus (AITC)

1. To induce scrunching using a chemical, e.g., the TRPA1 agonist AITC²⁸, planarians are ideally immersed in a bath of the chemical. If necessary, pipetting can be applied as described in section 2.1.2.3.

CAUTION: AITC is flammable, acutely toxic, can cause skin and eye irritation, respiratory and skin sensitization, and is hazardous to aquatic life. AITC oil should be handled in a fume hood. Prior to making stock solutions of AITC, put on appropriate PPE (nitrile gloves and a lab coat) and set up appropriate solid and liquid hazardous waste disposal containers.

2. In a fume hood, make a 10 mM stock solution of AITC in planarian water in a 50 mL centrifuge tube. This stock solution is useable for up to one month when stored at 4°C .

1. From this stock, prepare a 25 mL working solution of 100 μ M AITC in planarian water in a 50 mL centrifuge tube. This 100 μ M AITC solution will be used to induce scrunching in planarians.

NOTE: 100 μ M AITC induces consistent scrunching in *D. japonica* and *S. mediterranea* planarians²⁸. For other aquatic planarians, 100 μ M can serve as a starting concentration and can be adjusted accordingly.

2. Set up the experimental setup (refer to Section 1.1). Fill the arena with the AITC working solution and place it in a secondary container. The secondary container should hold at least twice the volume of the arena.

NOTE: Experiments can be carried out inside a fume hood for extra safety.

3. Transfer up to 10 planarians to the center of the arena and begin recording.
4. Once the planarians become desensitized and cease scrunching, stop recording. Remove the planarians from the AITC solution and rinse (refer to Section 1.1). Dispose of solid and liquid AITC waste in appropriate waste containers.
5. Verify the specificity of the response to AITC using RNAi to TRPA1²⁸ following standard protocols.

Representative Results

Extraocular near-UV perception in *S. mediterranea* planarians is TRPA1-dependent and has been proposed to be linked to H₂O₂ release¹⁷. Because H₂O₂ exposure induces TRPA1-dependent scrunching in *S. mediterranea* and *D. japonica* planarians²⁸, the steps in Section 2.1.4 can be used to test whether near-UV light exposure induces scrunching in

both species. While *D. japonica* planarians scrunch (10/10) when exposed to near-UV light, *S. mediterranea* planarians either exhibit tail thinning (7/10) as previously described¹⁷ or no response (3/10) (**Figure 4A,4B**). A quantification of the scrunching parameters, as outlined in Section 1.4, for the *D. japonica* planarians that exhibited at least 3 consecutive straight-line scrunches reveals characteristic scrunching parameters for this species^{7, 28} ($v_m = 0.84 \pm 0.14$, $|\Delta\epsilon|_{max} = 0.56 \pm 0.06$, $v^*_m = 0.47 \pm 0.07$, and $f_{elong} = 0.56 \pm 0.03$, values reported as mean \pm standard deviation for N=7).

In contrast, exposure to 250 μ M cinnamaldehyde, a known TRPA1 agonist in mice³⁴, causes scrunching in *S. mediterranea*^{7, 28} ($v_m = 0.46 \pm 0.08$, $|\Delta\epsilon|_{max} = 0.36 \pm 0.08$, $v^*_m = 0.16 \pm 0.04$, and $f_{elong} = 0.58 \pm 0.04$, values reported as mean \pm standard deviation for N=8) (**Figure 5A**), whereas *D. japonica* planarians at the same (and 1.6x the concentration) display a mixture of snake-like and oscillatory motion, interrupted by gliding and/or vigorous head turns (**Figure 5A**). A quantification of the (8/24) samples with at least three consecutive oscillations yields significantly lower values for 3 out of 4 parameters than expected for scrunching in this species ($v_m = 0.43 \pm 0.08$, $|\Delta\epsilon|_{max} = 0.39 \pm 0.03$, $v^*_m = 0.17 \pm 0.02$, and $f_{elong} = 0.54 \pm 0.06$, values reported as mean \pm standard deviation for N=8). Thus, while *D. japonica* appear to scrunch upon cinnamaldehyde exposure, a comparison of the calculated parameters with the literature values for this species^{7, 28} shows that the observed oscillatory motion is not scrunching. This example highlights the importance of quantitative measurements in conjunction with careful inspection of the raw behavioral data to properly interpret observed behaviors.

RNAi confirms the specificity of scrunching in response to cinnamaldehyde exposure in *S. mediterranea*. Within 180 seconds of exposure to 250 μ M cinnamaldehyde in planarian water 15/15 *unc22* (control) RNAi *S. mediterranea* planarians scrunched, whereas 0/16 *SmTRPA1* RNAi

planarians scrunched (**Figure 5B**), demonstrating that *S. mediterranea* scrunching in cinnamaldehyde requires *SmTRPA1*. Knockdown of *SmTRPA1* was confirmed through a 60 second exposure to a 100 μ M AITC bath²⁸.




Figure 1: Planarian behavior experimental setup.

(A) Sample experimental setup for studying planarian behavior. (B) 100 mm Petri dish arena centered in the field of view of the camera. [Please click here to view a larger version of this figure.](#)




Figure 2: Representative examples of the Fiji image analysis of planarians in arena.

(A) Selected region of interest, encompassing the full planarian path, indicated by the yellow rectangle. (B) Sample frames from the region of interest after duplication. (C) Subtracting the planarian from background and noise via thresholding (i) 8-bit image of planarian with noise, denoted by the asterisk. (ii) Binarized image of planarian after thresholding. (iii) Mask of planarian after setting filtering by size to remove noise. [Please click here to view a larger version of this figure.](#)




Figure 3: Plotting planarian length with respect to time.

(A) Raw plot of planarian length versus time for a scrunching *S. mediterranea* planarian. The asterisk denotes a moment when the planarian turned while scrunching. (B) Possible ways to trim scrunching data. (i) A correctly trimmed plot that removes the turning event data. (ii) An incorrectly trimmed plot that does not remove the turning event data. [Please click here to view a larger version of this figure.](#)




Figure 4: Species specific responses to near-UV light.

(A) Sample frames of *D. japonica* scrunching and *S. mediterranea* tail thinning in response to near-UV light. (B)

Representative oscillation plots of *S. mediterranea* and *D. japonica* in response to near-UV light. [Please click here to view a larger version of this figure.](#)




Figure 5: Species specific response to 250 μ M cinnamaldehyde, a TRPA1 agonist.

(A) Representative oscillation plots for *D. japonica* and *S. mediterranea* planarians in a 250 μ M cinnamaldehyde bath. (B)

Representative oscillation plots showing loss of scrunching in 250 μ M cinnamaldehyde in *SmTRPA1* RNAi *S. mediterranea* planarians. [Please click here to view a larger version of this figure.](#)

Supplemental Materials. [Please click here to download these materials.](#)

Discussion

Using this protocol, one can quantitatively study the effects of physical and chemical stimuli^{7, 28, 29} or genetic manipulation (RNAi)^{28, 29} on planarian locomotion. To maximize spatial

resolution, it is best to move the camera as close as possible to the arena while ensuring the entire arena is in the field of view. To increase throughput, the behavior of multiple planarians can be screened at once by recording multiple planarians simultaneously. When screening more than one planarian in a single arena, regions of interest can be drawn in Fiji to isolate individual planarians as described here or more advanced multi-object tracking can be employed. One issue with having multiple planarians in the same arena is that they can cross paths. This problem can be solved through the use of multi-well plates to isolate planarians from each other while still enabling simultaneous recording of many individuals to quantify behavior^{23, 24}. However, planarians will spend relatively more time at the wall in smaller arenas, requiring adjustments to the image analysis and limiting the resolution for scrunching/peristalsis quantification.

When stimuli are administered locally (e.g., pipetting⁷, amputation^{7, 28}, laser pointer¹⁷), it is crucial that the planarians are consistently stimulated in the same region because stimulating other body regions can potentially induce different behaviors. Different methods of delivery (such as pipetting or bath of a chemical) can also affect the consistency of the behavioral phenotype. Additionally, planarians can desensitize quickly²⁸, which needs to be taken into consideration when planning experiments as the same planarians should not be immediately reused for multiple experiments, either using the same or different stimuli. Finally, as shown here for near-UV exposure and cinnamaldehyde, it is important to be aware that the same stimulus can induce distinct behaviors in different planarian species. *D. japonica* scrunched when stimulated with near-UV light near the tail tip, while *S. mediterranea* planarians displayed tail thinning. In contrast, cinnamaldehyde exposure induced scrunching in *S. mediterranea* but not in *D.*

japonica planarians. Thus, while scrunching is a conserved response of various planarian species to noxious stimuli⁷, it has species specific parameters^{7, 28}, sensitivities²⁸, and inducers²⁸. Therefore, for a new species for which scrunching has not yet been parameterized, it is best to start with a well-conserved inducer, such as amputation⁷, to determine the species-specific parameters before testing the response to other stimuli.

One limitation of the analysis described here is that it does not account for turns and/or mixed behaviors, such as intermittent scrunching with head wiggling, gliding, or other body shape changes. However, close inspection of the raw data can help mitigate these issues if these instances are manually excluded from the analysis, as demonstrated in **Figure 3**. In addition, it is possible to add body shape analysis to the center of mass and length tracking described here and expand the protocol to quantify these other planarian behaviors. Given that the analysis does not make any assumptions about the studied organism, the protocol could in principle also be applied to other organisms that show similar types of behaviors.

The method of quantifying the different planarian gaits and distinguishing scrunching from peristalsis, as described here, assumes no prior training in computational image analysis or behavioral studies and does not require specialized equipment or software. To facilitate protocol adaptation, example data is provided in the Supplemental Material. The ease of obtaining and culturing planarians, as well as the ability to record behaviors without specialized equipment, makes planarian behavioral studies broadly accessible to research across all levels, from primary school classrooms to academic labs. A modified version of this protocol has been successfully used in a teaching laboratory setting that was

primarily composed of freshmen and sophomore students and included both prospective STEM and non-STEM majors.

The combination of molecular (RNAi) and chemical tools with quantitative behavioral analysis, as described in this protocol, allow researchers to gain mechanistic insights into the molecular control of behavior. Such work has uncovered some of the key mediators and neuronal circuits involved in planarian gliding^{19,20}, phototaxis^{17,35,36}, thermotaxis^{9,37}, and scrunching^{9,28,29}. Although planarian behaviors may not have direct corollary behaviors in higher organisms, such as humans, these behaviors represent fundamental neuronal functions important to all organisms - the ability to sense and process specific stimuli and react appropriately. Because of the conservation of key neuronal functions across different organisms, mechanistic studies in planarians can teach us more broadly about neuronal control of behavior. Additionally, analyzing planarian behavior in response to chemical exposure can be used to study the chemical's effects on the planarian nervous system^{23,24,25}, which may inform on potential risks to the human brain. In particular, scrunching induced by noxious heat was found to be a sensitive and specific endpoint for assaying neurotoxicity, because it becomes disrupted by exposure to certain classes of chemicals^{22,24,25,30}. Finally, the planarian's unique regenerative capabilities allow researchers to dissect the dynamics of how different behaviors are restored during neuroregeneration.

Disclosures

The authors have nothing to disclose.

Acknowledgments

The authors thank Mr. Tapan Goel for comments on the manuscript. This work was funded by NSF CAREER Grant 1555109.

References

1. Rink, J. C. Stem cell systems and regeneration in planaria. *Development Genes and Evolution*. **223**, 67–84 (2013).
2. Reddien, P. W., Alvarado, A. S. Fundamentals of Planarian Regeneration. *Annual Review of Cell and Developmental Biology*. **20**, 725–757 (2004).
3. Cebrià, F. Regenerating the central nervous system: how easy for planarians! *Development Genes and Evolution*. **217**, 733–748 (2007).
4. Pearl, R. Memoirs: The Movements and Reactions of Fresh-Water Planarians: A Study in Animal Behaviour. *Journal of Cell Science*. **s2-46** (1903).
5. Mc Connell, J. *A Manual of Psychological Experimentation on Planarians*. Planarian Press, (1967).
6. Talbot, J., Schötz, E.M. Quantitative characterization of planarian wild-type behavior as a platform for screening locomotion phenotypes. *Journal of Experimental Biology*. **214**, 1063–1067 (2011).
7. Cochet-Escartin, O., Mickolajczk, K. J., Collins, E.M. S. Scrunching: a novel escape gait in planarians. *Physical Biology*. **12**, 055001 (2015).
8. Inoue, T. et al. Planarian shows decision-making behavior in response to multiple stimuli by integrative brain function. *Zoological Letters*. **1**, 1–15 (2015).
9. Arenas, O. M. et al. Activation of planarian TRPA1 by reactive oxygen species reveals a conserved mechanism

- p>for animal nociception.
- Nature Neuroscience*
- .
- 20**
- , 1686–1693 (2017).
10. Shomrat, T., Levin, M. An automated training paradigm reveals long-term memory in planarians and its persistence through head regeneration. *Journal of Experimental Biology*. **216**, 3799 LP – 3810 (2013).
11. Blackiston, D., Shomrat, T., Nicolas, C. L., Granata, C., Levin, M. A Second-Generation device for automated training and quantitative behavior analyses of Molecularly-Tractable model organisms. *PLoS One*. **5**, 1–20 (2010).
12. Ross, K. G., Currie, K. W., Pearson, B. J., Zayas, R. M. Nervous system development and regeneration in freshwater planarians. *Wiley Interdisciplinary Reviews-Developmental Biology*. **6**, e266 (2017).
13. Cebrià, F. et al. The expression of neural-specific genes reveals the structural and molecular complexity of the planarian central nervous system. *Mechanisms of Development*. **116**, 199–204 (2002).
14. Mineta, K. et al. Origin and evolutionary process of the CNS elucidated by comparative genomics analysis of planarian ESTs. *Proceedings of the National Academy of Sciences of the United States of America*. **100**, 7666–71 (2003).
15. Ross, K. G. et al. SoxB1 Activity Regulates Sensory Neuron Regeneration, Maintenance, and Function in Planarians. *Developmental Cell*. **47**, 331–347.e5 (2018).
16. Nishimura, K. et al. Reconstruction of Dopaminergic Neural Network and Locomotion Function in Planarian Regenerates. *Developmental Neurobiology*. **67**, 1059–1078 (2007).
17. Birkholz, T. R., Beane, W. S. The planarian TRPA1 homolog mediates extraocular behavioral responses to near-ultraviolet light. *Journal of Experimental Biology*. **220**, 2616–2625 (2017).
18. Currie, K. W., Molinaro, A. M., Pearson, B. J. Neuronal sources of hedgehog modulate neurogenesis in the adult planarian brain. *Elife*. **5** (2016).
19. Talbot, J. A., Currie, K. W., Pearson, B. J., Collins, E.M. S. Smed-dynA-1 is a planarian nervous system specific dynamin 1 homolog required for normal locomotion. *Biology Open*. 1–8 (2014).
20. Currie, K. W., Pearson, B. J. Transcription factors *lhx1/5-1* and *pitx* are required for the maintenance and regeneration of serotonergic neurons in planarians. *Development*. **140**, 3577–88 (2013).
21. Hagstrom, D. et al. Planarian cholinesterase: molecular and functional characterization of an evolutionarily ancient enzyme to study organophosphorus pesticide toxicity. *Archives of Toxicology*. **92**, 1161–1176 (2018).
22. Hagstrom, D., Cochet-Escartin, O., Collins, E.M. S. Planarian brain regeneration as a model system for developmental neurotoxicology. *Regeneration*. **3**, 65–77 (2016).
23. Hagstrom, D., Cochet-Escartin, O., Zhang, S., Khuu, C., Collins, E.M. S. Freshwater planarians as an alternative animal model for neurotoxicology. *Toxicological Sciences*. **147**, 270–285 (2015).
24. Zhang, S., Hagstrom, D., Hayes, P., Graham, A., Collins, E.M. S. Multi-behavioral endpoint testing of an 87-chemical compound library in freshwater planarians. *Toxicological Sciences*. 26–44 (2019).

25. Zhang, S., Hagstrom, D., Siper, N., Behl, M., Collins, E.M. S. Screening for neurotoxic potential of 15 flame retardants using freshwater planarians. *Neurotoxicology and Teratology*. **73**, 54–66 (2019).
26. Wu, J. P., Li, M. H. The use of freshwater planarians in environmental toxicology studies: Advantages and potential. *Ecotoxicology and Environmental Safety*. **161**, 45–56 (2018).
27. Rompolas, P., Azimzadeh, J., Marshall, W. F., King, S. M. Analysis of ciliary assembly and function in planaria. *Methods in Enzymology*. **525**, 245–264 (2013).
28. Sabry, Z. et al. Pharmacological or genetic targeting of Transient Receptor Potential (TRP) channels can disrupt the planarian escape response. *PLoS One*. 753244 (2019).
29. Cochet-Escartin, O., Carter, J. A., Chakraverti-Wuerthwein, M., Sinha, J., Collins, E. M. S. Slo1 regulates ethanol-induced scrunching in freshwater planarians. *Physical Biology*. **13**, 1–12 (2016).
30. Hagstrom, D., Truong, L., Zhang, S., Tanguay, R. L., Collins, E.M. S. E.M. S. E.M. S. Comparative analysis of zebrafish and planarian model systems for developmental neurotoxicity screens using an 87-compound library. *Toxicological Sciences*. kfy180 (2019).
31. Schindelin, J. et al. Fiji: an open-source platform for biological-image analysis. *Nature Methods*. **9**, 676–682 (2012).
32. Akiyama, Y., Agata, K., Inoue, T. Spontaneous Behaviors and Wall-Curvature Lead to Apparent Wall Preference in Planarian. *PLoS One*. **10**, e0142214 (2015).
33. Paskin, T. R., Jellies, J., Bacher, J., Beane, W. S. Planarian Phototactic Assay Reveals Differential Behavioral Responses Based on Wavelength. *PLoS One*. **9**, e114708 (2014).
34. Petrus, M. et al. A role of TRPA1 in mechanical hyperalgesia is revealed by pharmacological inhibition. *Molecular Pain*. **3**, 40 (2007).
35. Takano, T. et al. Regeneration-dependent conditional gene knockdown (Readyknock) in planarian: Demonstration of requirement for Djsnap-25 expression in the brain for negative phototactic behavior. *Development, Growth & Differentiation*. **49**, 383–394 (2007).
36. Nishimura, K. et al. Identification of glutamic acid decarboxylase gene and distribution of GABAergic nervous system in the planarian *Dugesia japonica*. *Neuroscience*. **153**, 1103–14 (2008).
37. Inoue, T., Yamashita, T., Agata, K. Thermosensory signaling by TRPM is processed by brain serotonergic neurons to produce planarian thermotaxis. *Journal of Neuroscience*. **34**, 15701–14 (2014).